



Configuring DTSA-II

Nicholas W. M. Ritchie
Physicist, Microanalysis Research Group
National Institute of Standards and Technology
Gaithersburg, MD 20899-8371
nicholas.ritchie@nist.gov








Spectrum Simulation
via Monte Carlo or
analytical models

When DTSA-II simulates
spectra, it uses the
same instrument,
detector and calibration
definitions as when you
read in spectra or when
you quantify spectra.

Not all of the meta-data in
the detector definitions is
important for all uses.
Simulation tends to use more
of the meta-data items than
quantification. However,
there is one shared definition
for both making easy to
“simulate what you measure”.







But I don't know some of the values you ask about my detector...

Don't fret. It is likely the vendor who sold you your detector doesn't have good answers for every question. Even if they did, they don't feel compelled to tell us. Some claim that it would put them at a competitive disadvantage if they told us. More really don't know. The vendor that they purchase their detector crystals won't tell them.


Fortunately, for almost every purpose, it really doesn't matter.

(This seems to bother some OC people.)

So why do I continue to ask? Maybe your purpose is special and it does matter. If this is the case (and you will know), then DTSA-II will allow you to specify the details.

Personally, I don't sweat it too much. Usually, I take most of the default values unless I know otherwise. Some parameters do matter and...

I will point out those parameters that do matter.



NIST DTSA-II

File Process Tools Report Help

Counts

Same scale
899 x 259

Cr2N_GLR_1_10kV20nA_Med_54kHz6DT_100s
Cr2N_GLR_2_10kV20nA_Med_54kHz6DT_100s
Cr2N_GLR_3_10kV20nA_Med_54kHz6DT_100s

NIST DTSA-II - Preferences

User Information

Give me a name!

Give the instrument a user friendly name.

Instrument name Mega SEM ++

Min beam energy (keV) 5.0 Max beam energy (keV) 30.0

Add Si(Li) detector Add SDD detector
Add microcalorimeter Import detector

Then add a detector by creating a conventional Si(Li) detector, a modern Silicon Drift Detector (SDD) or by importing an detector definition exported by someone else.

Probe
Oxford
FEI FIB
Personal SEM
842.10 JEOL
Give me a name!

Default Detector
JXA-8500F
Bruker 5 eV/ch

Spectrum List
Cr2N_GLR_1_10kV20nA
Cr2N_GLR_2_10kV20nA
Cr2N_GLR_3_10kV20nA

OK Cancel Apply

Instrument
JXA-8500F

3 spectra selected

Atomic Fraction

NIST DTSA-II

File Process Tools Report Help

Counts

Same scale
891 x 259

60 000
40 000
20 000
0

5 10
15 20

Spectrum Report

Default Detector
JXA-8500F
Bruker 5 eV/ch

NIST DTSA-II - Preferences

- User Information
- Quantitative algorithms
- Instruments and Detectors
 - McCrone
 - Crap
 - FPS
 - MIRA-3
 - JXA-8500F
 - Probe
 - Oxford
 - FEI FIB
 - Personal SEM
 - 842.10 JEOL
 - Mega SEM ++
 - Detector - HyperSDD (medium TC)**

Detector - HyperSDD (medium TC)

Status

Enable detector ☒

Name

Detector name HyperSDD (medium TC, 2048 ch)

Globally Unique Identifier f354c815-7d62-daa1-891a-4d473a68929a

Import

Import from spectrum

Often vendor's spectrum files will contain some of the data items in this panel. Use import to populate the panel with the vendor's values. Use care. Vendor's values are sometimes wrong.

Scroll to see more parameters.

Atomic Fraction

3 spectra selected

A new detector is added to your instrument.

Give the detector a descriptive name.

You are presented with a panel to edit the detector definition.

You are likely to be able to populate many of the detector parameters from a spectrum data file collected on this detector. Be careful however. Some of the data in a vendors spectrum file may be incorrect. Imported meta-data is identified by a yellow background.

NIST DTSA-II

File Process Tools Report Help

60 000 Same scale 89% x 259

Cr2N_GLR_1_10kV20nA_Med_54kHz6DT_100s
Cr2N_GLR_2_10kV20nA_Med_54kHz6DT_100s
Cr2N_GLR_3_10kV20nA_Med_54kHz6DT_100s

Specify the vacuum tight window that protects the detector when the SEM chamber is vented. (Critical)

5 10 15 20

Spectrum Report

Default Detector JXA-8500F

Bruker 5 eV/ch

Spectrum List

- Cr2N_GLR_1_10kV20nA
- Cr2N_GLR_2_10kV20nA
- Cr2N_GLR_3_10kV20nA

JXA-8500F

- Probe
- Oxford
- FEI FIB
- Personal SEM
- 842.10 JEOL
- Mega SEM ++
- Detector - HyperSDD (medium TC, ...)

Detector - HyperSDD (medium TC, ...)

Window

Window

Position

Side View

Electron Column

Detector snout

Sample

Working Distance

Elevation

Elevation angle 35.0

Window

- No window
- Beryllium (25 μm)
- Moxtek AP 1.3
- Moxtek AP 1.7
- Moxtek AP 3.3 (model)
- Moxtek AP 3.3 (manufacturer's table)
- No window
- Diamond (0.45 μm)
- Boron Nitride (0.25 μm)

Electron Column

Azimuth

Edit the properties of this detector.

Atomic Fraction

3 spectra selected

The window is the single most important parameter for simulation as it has the strongest influence on the low energy efficiency of the detector.

Fortunately, almost all modern detectors have a "light element" polymer window similar to the AP 3.3. If you can see x-rays from carbon and oxygen, you have a "light element" window.

NIST DTSA-II

File Process Tools Report Help

60 000 Same scale 891 x 259

DTSA-II provides visual guidance to help you to interpret the various different geometric parameters.

What DTSA-II correctly refers to as "elevation" is often called "take-off angle." Actually, "take off angle" is elevation plus a component of the sample tilt modified by working distance.

Detector

Position

Side View

Top View looking down

Electron Column

Detector snout

Sample

Working Distance

Elevation

Electron Column

Azimuth

Elevation angle 35.0

Azimuthal angle 0.0

Optimal working distance 20.0

Edit the properties of this detector.

Elevation angle is the most critical parameter for quantification. (Critical)

Azimuthal angle is important when modeling shaped samples. You need to know where the detector is to interpret the emission images. (Important for simulation)

3 spectra selected

Cr2N_GLR_1_10kV20nA

Cr2N_GLR_2_10kV20nA

Cr2N_GLR_3_10kV20nA

Default Detector

JXA-8500F

Brucker 5 eV/ch

Spectrum List

JXA-8500F

Probe

Oxford

FEI FIB

Personal SEM

842.10 JEOL

Mega SEM ++

Detector - HyperSDD (medium TC)

-481.9 eV

5.01 eV/channel

0 nm

JXA-8500F

NIST DTSA-II

File Process Tools Report Help

60 000 Same scale 891 x 259

Cr2N_GLR_1_10kV20nA_Med_54kHz6DT_100s
Cr2N_GLR_2_10kV20nA_Med_54kHz6DT_100s
Cr2N_GLR_3_10kV20nA_Med_54kHz6DT_100s

NIST DTSA-II - Preferences

Detector - HyperSDD (medium TC)

Elevation angle 35.0 °
Azimuthal angle 0.0 °
Optimal working distance 20.0 mm
Sample-to-detector distance 50.0 mm

Crystal parameters

Detector Area 10.0 mm²
Gold layer 0.0 nm
Aluminum layer 10.0 nm
Nickel layer 0.0 nm
Dead layer 0.000 nm

Edit the properties of this detector.

OK Cancel

3 spectra selected

The “sample-to-detector distance” is the distance from the intersection of cylindrical axis of the detector with the electron beam to the front face of the detector crystal.
You may take a guess and adjust this value until the intensity of simulated spectra match measured spectra.
(Important)

The “optimal working distance” is the distance, as measured by the objective lens, at which the cylindrical axis of the x-ray detector intersects the electron beam.
(Important)

Your vendor should provide the “detector area.” (Critical)

These represent various different conductive coatings layered on the front of the detector crystal to make it conductive.
(Less important)

The “dead layer” is an inactive layer of silicon on the front face of the x-ray detector. It absorbs some x-rays but does not produce photo-electrons to measure.
(Less important)

NIST DTSA-II

File Process Tools Report Help

60 000 Same scale 89 x 259

NIST DTSA-II - Preferences

Detector - HyperSDD

Dead layer 0.000 μm

Thickness 0.45 mm

Configuration

Number of channels 2048 channels

Zero strobe discriminator 0.0 eV

Base Performance

Energy scale 5.0 eV/channel

Zero offset 0.0 eV

Resolution 128.0 eV FWHM at Mn K α

Export detector

Edit the properties of this detector.

OK

3 spectra selected

The “thickness” of the detector crystal influences the efficiency of the detector at photon energies above 10 keV. (Less important)

The number of channels of spectrum data recorded in a spectrum file. (Important)

To a good approximation, the mapping between detector channel and x-ray energy is linear. The scale and offset allow you to specify the slope and offset of the mapping. Enter nominal values here. The exact values can be calibrated. (Critical)

The resolution is a detector performance metric. It is the full width half maximum of a measured Mn K α peak. Enter a nominal value here. The exact value can be calibrated. (Important)


Some vendors inject a fake peak at 0 eV as a mechanism 1) to stabilize the zero offset; 2) to measure the amplifier noise characteristics. The “zero strobe discriminator” allows you to specify a low energy below which all counts are assumed to be associated with the zero strobe. (Important)

This button will export a full detector definition to a file. This allows users to import other people’s detector definitions from the “Instrument” page.

Parameters by importance

- Critical
 - “Window”, “Elevation angle”, “Detector area”, “Energy scale”, “Zero offset”
- Important
 - “Azimuthal angle”, “Optimal working distance”, “Sample-to-detector distance”, “Number of channels”, “Zero strobe discriminator”, “Resolution”
- Less important
 - “Gold layer”, “Nickel layer”, “Aluminum layer”, “Dead layer”, “Thickness”

Do your best. It will
probably be good enough.



NIST DTSA-II

File Process Tools Report Help

Counts

Same scale
89% x 259

60 000
40 000
20 000
0

5 10
15 20

Spectrum Report Command

Default Detector
Mega SEM ++
HyperSDD (medium TC, 2048 ch)

Spectrum List
Cr2N_GLR_1_10kV20nA_Med_54kHz
Cr2N_GLR_2_10kV20nA_Med_54kHz
Cr2N_GLR_3_10kV20nA_Med_54kHz

Calibrate an EDS detector

Measured spectrum

Specify a spectrum
Spectrum ...

Specify the material
Material Mn standard
Live time sec. Probe current nA
Fit type Linear (default)

Specify an effective date
Effective date

Message: More...

Back Next Finish Cancel

When you create a detector, you enter a nominal resolution, scale and zero offset. The “calibration alien” allows you refine this estimate using a measured spectrum.

The calibration process is described in another presentation.

3 spectra selected

Model	Moxtek AP 3.3 (manufacturer's table)
Angle	40°
Offset	-481.9 eV
Resolution	5.01 eV/channel
Thickness	0 nm
Detector	JXA-8500F

Mass Fraction Atomic Fraction