

Quantifying STEM-mode spectra using ζ -factors


Getting ready

You will need the following information to perform a STEM-mode quantification. STEM, of course, means scanning transmission electron microscopy and the samples we are talking about are unsupported thin-films or particles on an essentially transparent thin film – no substrate.

1. Standard spectra collected from thin films of known composition and mass-thickness for each element in the unknown.
 - a. Maybe the single largest challenge in STEM quant is finding suitable artifacts to act as standards. Usually, standard thin-films are prepared and then characterized by other techniques.
 - b. Standards may be single element or multi-element but there may not be any interferences between elements in the standard.
 - c. The beam dose (live time · probe current) must be known for each standard spectrum
 - d. The standard thin-films should be thin enough that less than a percent or two of the incident beam energy is lost in the film.
 - e. The standards should have been collected at the same beam energy as the unknown spectrum,
2. Reference spectra collected from bulk or thin samples for any spurious peaks like Cu from the grid.
 - a. The mass-thickness or the probe dose is not necessary for reference spectra.
3. The unknown spectrum along with the probe dose.
 - a. The unknown film should be thin enough that less a percent or two of the incident beam energy is lost in the film.
 - b. If the unknown and standard spectra are collected without interruption on a stable instrument, it may not be necessary to measure the probe current so long as it remains constant.

Prepare the standard spectra

1. It is easiest if standard spectra contain all the information required to facilitate the quantification. DTSA-II provides a mechanism to ensure that this data is encapsulated within the DTSA-II-extended EMSA spectrum files. Your spectra files should contain data specifying the:
 - a. Probe current
 - b. Live time
 - c. Composition
 - d. Mass-thickness
 - e. Beam energy
 - f. Detector
2. Fortunately DTSA-II provides the “make ‘Standard Bundle’” dialog to ensure that spectra contain all the necessary data.




The “make standard” dialog is accessed by selecting a spectrum (or multiple spectra if you wish to combine them) in the “spectrum list” and selecting the “make standard” option from the “spectrum list” right-click popup menu. Since we wish to perform thin film quant, we will need to specify a material with a density and the thickness of the film.

3. Use the “save” button in the lower, left of the “make standard” dialog to write the spectrum to disk. Remember where you saved it as you will need to reload this spectrum from within the “quantification wizard.”
4. Both the unknown and standards should be associated with the same instrument and detector.

Prepare the unknown spectra

1. It is also helpful to prepare the unknown spectra. The data required for the unknown include:
 - a. Probe current
 - b. Live time
 - c. Beam energy
 - d. Detector
2. You can enter this information using the “edit spectrum properties” dialog. This is accessed from the “spectrum list” right-click popup menu or from the “Tools” item on the main menu.

Again you should select the spectrum you wish to edit in the “spectrum list.”



3. Ensure that the data on the “conditions” tab is correct and select “Ok.” You can write the spectrum to disk as an EMSA file but you should also retain the unknown spectrum in the “spectrum list.”


Preparing reference spectra

1. Reference spectra are used to strip spurious elements which are common in STEM spectra. Since the spectra are only used for the peak shape, bulk or thin film spectra will work and very little information is actually required. The only additional piece of information that it is useful to include in a reference spectrum is the composition of the material from which it was collected.
2. The “assign material” menu item in the “spectrum list” popup menu can be used to assign a material to the reference. This helps the program know which elements the spectrum can be used as a reference for.

Quantifying the unknown


1. Select one or more edited thin-film spectra to quantify in the “spectrum list”.

2. Select the “quantification alien...” menu item from the “tools” main menu.



3. Select the “quantify a STEM spectrum using MLLSQ fitting and ζ -factors” option.

4. Ensure the instrument, detector, calibration and beam energy are correct. They should be as this information is read from the s



The image shows a software window titled "Quantification Alien" with a close button (X) in the top right corner. Inside the window, there is a header section with a grey alien head icon on the left. To the right of the icon, the text "Specify the instrument" is displayed in a large, bold font. Above this text, it says "Previous: Select a quantification mode", and below it, "Next: Specify standard spectra (STEM)".

The main area of the window contains four labeled sections, each with a dropdown menu:

- Instrument**: A label followed by a horizontal line, then "Acquired on the" followed by a dropdown menu showing "MIRA-3".
- Detector**: A label followed by a horizontal line, then "using the" followed by a dropdown menu showing "SDD (Medium, 4096)".
- with calibration**: A label followed by a horizontal line, then a dropdown menu showing "FWHM[Mn K α]=128.0 eV - initial".
- Setting**: A label followed by a horizontal line, then "at a beam energy of" followed by a text input field containing "30.0" and the unit "keV".

At the bottom of the window, there is a "Message:" label on the left and a "More..." button on the right. Below these, there are four buttons: "Back", "Next", "Finish", and "Cancel". The "Cancel" button is highlighted with a blue border.

5. The next panel allows you to assign standard spectra to elements. Use the "add" and "remove" buttons to open standard spectra from disk. The "edit" button allows you to edit the properties

of the selected standard spectrum.

Quantification Alien

Previous: *Specify the instrument*

Specify standard spectra (STEM)

Next: *Other elements to strip*


Spectrum	Composition	Elements	Dose	Mass-Thickness
Ar std	Ar	Ar	100.0 nA·s	12.5 $\mu\text{g}/\text{cm}^2$
Fe std	Fe	Fe	200.0 nA·s	23.6 $\mu\text{g}/\text{cm}^2$
MgO std	MgO	O, Mg	60.0 nA·s	35.8 $\mu\text{g}/\text{cm}^2$
Si std	Si	Si	80.0 nA·s	16.3 $\mu\text{g}/\text{cm}^2$
CaF2 std	CaF2	Ca	120.0 nA·s	15.9 $\mu\text{g}/\text{cm}^2$

Add
Remove
Edit

Message: More...

Back Next Finish Cancel

6. You may be asked to specify which element to associate multi-element spectra with. In this case, MgO could be a standard for either Mg or O or both Mg and O. Type in the elements separated by commas or spaces in the “elements” edit box or use the “...” button to select the elements from a periodic table.



Typing "Mg, O" assigns the "MgO std" to both elements.

- The next panel allows you to assign spectra to strip zero or more elements from the spectrum. Stripping elements is particularly important if a spurious characteristic peak from a mounting grid or a sample support is present in the spectrum and this peak interferes with an element you wish to measure.

Quantification Alien

Previous: *Specify standard spectra (STEM)*

Next: *Specify the unknown spectra*

Other elements to strip

Element	Spectrum
Titanium	Ti std
Copper	Cu std

Add

Remove

Message: Specify the elements to strip

More...

Back Next Finish Cancel

My STEM has Ti in the holder and Cu in the grid which are visible in the measured spectrum.

8. The next panel allows you to review the spectra being quantified and edit the spectrum properties.

Quantification Alien

Previous: *Other elements to strip*

Specify the unknown spectra

Next: *Quantification Results*

Spectrum	Dose	Beam Energy
SRM-2063a std with scatter	240.2 nA's	30.0 keV

Remove


Edit

Message: Specify the unknown spectra

More...


Back Next Finish **Cancel**

9. The final panel performs the quantification and summarizes the results.



Select the "finish" button to exit the "quantification alien."

10. In addition to the original spectra, you will find a second spectrum corresponding to the residual spectrum in the "spectrum list." Examining the residual spectrum is a good way to determine if the quantification account for all the elements in the unknown.



This residual shows that we've successfully accounted for all the elements in the unknown including the unlabeled Ti and Cu instrument scatter peaks.

11. The “report” tab contains detailed description of the quantification and the results.

Conditions

Item	Value
Instrument	MIRA-3
Detector	SDD (Medium, 4096)
Beam Energy	30.0 keV
Correction Algorithm	ζ -factors
Mass Absorption Coefficient	NIST-Chantler 2005
Stripped elements	Ti, Cu

The conditions table summarizes many instrumental parameters.

Standards

Element	Material				Dose	Mass-Thickness $\mu\text{g}/\text{cm}^2$	Spectrum
O	MgO				60.000	35.800	MgO std
	Element	Mass Fraction	Mass Fraction (normalized)	Atomic Fraction			
	O	0.3970	0.3970	0.5000			
	Mg	0.6030	0.6030	0.5000			
Mg	MgO				60.000	35.800	MgO std
	Element	Mass Fraction	Mass Fraction (normalized)	Atomic Fraction			
	O	0.3970	0.3970	0.5000			
	Mg	0.6030	0.6030	0.5000			
Si	Si				80.000	16.310	Si std
	Element	Mass Fraction	Mass Fraction (normalized)	Atomic Fraction			
	Si	1.0000	1.0000	1.0000			
Ar	Ar				100.000	12.500	Ar std
	Element	Mass Fraction	Mass Fraction (normalized)	Atomic Fraction			
	Ar	1.0000	1.0000	1.0000			
Ca	CaF2				120.000	15.900	CaF2 std
	Element	Mass Fraction	Mass Fraction (normalized)	Atomic Fraction			
	F	0.4867	0.4867	0.6667			
	Ca	0.5133	0.5133	0.3333			
Fe	Fe				200.000	23.610	Fe std
	Element	Mass Fraction	Mass Fraction (normalized)	Atomic Fraction			
	Fe	1.0000	1.0000	1.0000			

The standard panel summarizes the standard spectra and the materials from which they were collected.

Results

Spectrum	Quantity	O		Mg		Si		Ar		Ca		Fe		Sum
SRM-2063a std with scatter	Line	O All		Mg All		Si All		Ar K-family		Ca K-family		Fe Ka		
	k-ratios	0.7081	± 0.0059	0.0871	± 0.0013	0.3530	± 0.0030	0.0063	± 0.0010	0.3309	± 0.0059	0.1093	± 0.0032	
Mass-thickness 23.3±0.3 µg/cm²	mass fraction	0.4384	± 0.0017(K) 0.0008(MAC) [0.0018]	0.0808	± 0.0001(K) 7.3·10 ⁻⁵ (MAC) [0.0002]	0.2507	± 0.0012(K) 0.0003(MAC) [0.0012]	0.0034	± 0.0005(K) 5.0·10 ⁻⁵ (MAC) [0.0005]	0.1161	± 0.0005(K) 0.0002(MAC) [0.0005]	0.1106	± 0.0017(K) 0.0002(MAC) [0.0017]	1.0000
I = 1.000 nA	norm(mass fraction)	0.4384	± 0.0018	0.0808	± 0.0002	0.2507	± 0.0012	0.0034	± 0.0005	0.1161	± 0.0005	0.1106	± 0.0017	-
LT = 240.2 s	atomic fraction	0.6142	± 0.0026	0.0746	± 0.0001	0.2000	± 0.0010	0.0019	± 0.0003	0.0649	± 0.0003	0.0444	± 0.0007	
Residual	C:\Users\ritchie.NIST\AppData\Local\NIST\WIST DTSA-II Reports\2015\November\6-Nov-2015\residual\436558380414768936.msa													
Notes	Uncertainties are 1 σ and labeled by source. (K: k-ratio, MAC: mass absorption coefficient, []: combined)													

The results table provides a detailed results report including the measured k-ratios, the measured mass-thickness, the measured mass-fraction which is also expressed in normalized and atom fraction forms. The blue links allow the original spectrum and residual to be reloaded from disk.

The report provides uncertainties which are based on count statistics and assumed uncertainties in the mass absorption coefficients.

Nicholas W. M. Ritchie
6-Nov-2015