



Quantification

A brief look at how quantification of x-ray spectra is performed

What material produces a spectrum that looks like this?





==

?




+

?




+

?





+

?




+

k-ratios



$$k = \frac{35,310 \pm 373}{384,807 \pm 782} = 0.0917 \pm 0.0010$$



k-ratios



0.1076




0.1792




0.2947



0.0921




||




+



+




+



+







Matrix correction




$$C_A = \frac{\{Z \cdot A \cdot F\}_{Unk,A}}{\{Z \cdot A \cdot F\}_{Std,A}} k_A$$

Matrix corrections turn k-ratios into compositions

Monte Carlo models




Matrix correction



Results

Spectrum	Quantity	O	Mg	Si	Ca	Fe	Sum
K-411 3	Line		Mg All	Si All	Ca K-family	Fe Ka	
Bulk	Z · A · F	- -	1.024 0.738 1.007	0.982 0.719 1.002	0.986 1.000 1.006	0.842 0.991 1.000	
	k-ratios	- -	0.1076 ± 0.0004	0.1792 ± 0.0004	0.2947 ± 0.0013	0.0921 ± 0.0007	
D-H = 20.123 keV	mass fraction	0.4222 ± 0.0002	0.0853 ± 0.0002	0.2535 ± 0.0002	0.1144 ± 0.0002	0.1103 ± 0.0002	0.9858
I = 12.600 nA	norm(mass fraction)	0.4283 ± 0.0003	0.0866 ± 0.0003	0.2571 ± 0.0004	0.1161 ± 0.0005	0.1119 ± 0.0007	-
LT = 300.0 s	atomic fraction	0.6031 ± 0.0005	0.0802 ± 0.0002	0.2062 ± 0.0001	0.0653 ± 0.0002	0.0451 ± 0.0003	




Element	Standard	Measured	Difference	% Difference
Si	25.38	25.35	-0.03	-0.1
Fe	11.21	11.03	-0.18	-1.6
O	42.37	42.22	-0.14	-0.3
Ca	11.06	11.44	0.39	3.5
Mg	8.85	8.53	-0.31	-3.5
5	98.86	98.58	0.55	5.3

Assumptions

- Homogeneous material
- Flat
- Polished
- Clean
- Mounted normal to the electron beam
- Conductive or conductively coated
- Standards / unknown measured under same conditions

What did I gloss over?


- How were the samples prepared?
- Why did I select the standards I selected?
- Why did I select the beam energy I selected?
- Are any references required?
- Which line (K, L or M) should I use to quantify?
- How many x-rays do I need to collect (count time / probe dose)?
- Where do I get standard materials?



Geller



SPI



Custom

Why is it that we are willing to spend \$50k+ on a detector but not an additional \$10k on accessories to make our detector more accurate?